

# Biner LC

Light curing type Cavity Liner

## Light curing time 40s

Biner LC is a light cure, fluoride releasing, radiopaque cavity liner and base material specially formulated for use with adhesives and composites and with conventional restorative materials.

Biner LC contains Hydroxy calcium phosphate in a urethane dimethacrylate resin. It chemically bonds to adhesive primers, composites and other resin based material, and it micromechanically bonds to dentin.

### INDICATION

- Protect pulp of the tooth
- Lining under all filling materials
- Protection when applying the total etch technique

### CONTRAINDICATIONS

- Biner LC should not be used if patients are known to be allergic to methacrylate resins

### INSTRUCTIONS

- Isolate prepared teeth and air dry gently
- Place a disposable tip securely onto Biner LC syringe
- Apply Biner LC to the dentin in increments. Remove Biner LC from enamel margin and light cure for 40 seconds, depending on the light curing unit.
- Apply bonding system to enamel. Light cure. Continue with layer of composite until cavity is filled.

### SIDE EFFECTS

- None known to date

### CAUTIONS

- Avoid contact of Biner LC with the skin/mucous membrane or eyes
- None polymerized Biner LC may be causes lightly, reversible irritation or lead to a general sensitization against methacrylates.
- Keep away from child.
- For dental use only
- Do not use after the date of expiration.
- Exposure to light may be lead to premature polymerization

### STORAGE

- Store at room temperature
- Close syringe cap immediately after use.
- Avoid from directly sun light.



## INSTRUCTION FOR USE

### SHELF LIFE

- 3 years

### COMPONENTS

- Biner LC 2g×2 syringes

### SYMBOLS



Use By



Batch Code



0120

CE-Mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC



Authorized representative in the EU



Attention, See instruction for use



Manufacturer

Made in Korea

Manufactured by

**META BIOMED CO.,LTD.**

634, Yeonje-ri, Gangseo-myeon  
Cheongwon-gun, Chungcheongbuk-do, Korea  
Tel: +82-43-216-0433  
Fax: +82-43-217-1988  
E-mail: info@meta-biomed.com  
Web site: http://www.meta-biomed.com

DM-85-2

## CHECKING SHEET

BINER LC

TITLE

CODE DM-85-2

SIZE 157x189mm

PAPER 60g 모조

MAIN COLOR INDIGO 2767C 별색 단도

SPECIAL

FUNCTION

PRINT COLOR



INDIGO  
2767C

DATE 2009. 07. 10.

DESIGNER 진희남

RECOGNITION 김상일

CONFIRMATION META BIOMED CO.,LTD.

DATE

NAME

SIGN

→ 2011.05.17 DM-85-2 주소변경 (홍선미)