


MetaFix Knotless PK™



Order Information

— Anchor

Model No	Description
MA-4521-P1A	Φ4.5 x 20.4mm, Black handle, #2 (White&Blue)
MA-5521-P1A	Φ5.5 x 21.0mm, White handle, #2 (White&Blue)

— Instrument

Model No	Description
MI-38293-A	3.8mm awl

270 Osongsaengmyeong1-ro, Osong-eup,
Heungdeok-gu, Cheongju-si, Chungbuk, Korea

Phone: +82 43 230 8836 Fax: + 82 43 230 8969
medical@mata-biomed.com www.meta-biomed.com


HO-25-1

MetaFix Knotless PK™

Feature

- Lateral row anchor
- Tap-in insertion type
- Optimized suture-bridge technique
- Inner screw suture fixation mechanism
- Adjustable anchor suture retention
- Strong bone pullout strength

PEEK-OPTIMA™ reinforced offers enhanced mechanical and physical properties compared to the base polymer, making it suitable for more demanding load-bearing applications that are currently beyond the capability of PEEK-OPTIMA™ natural.



3.8mm awl

Secure Suture Management

Torque indicator