



Alondra

Suction Irrigator with cautery



The Alondra

Suction Irrigator with cautery device has 2 types of handle design.

Ergonomic Gun type handle design with easy control access.

Ergonomic handle design

- Improved gun type handle for optimized ergonomics enables reduce surgeon wrist fatigue.
- Easy control flow with finger-tip access.

Energy Button



ooth and shorten operation time.

Trumpet Type Handle

- 360° Rotating shaft
- 8° Angled shaft



Electrode Features

Teflon coating prevents tissue from sticking to the tip when during surgery.

Multiple Electrode Option

Multiple design tip options allows surgeon to pick best type of tip for operation.



L-hook



Spatula



Curved dissector



J-hook

Order Information

Model Name	Handle Type	Electrode Type	Shaft Length
Alondra-L	Gun	L Hook	350mm
Alondra-J		J Hook	
Alondra-S		Spatula	
Alondra-D		Curved Dissector	
Alondra-SL	Trumpet Straight	L Hook	
Alondra-SJ		J Hook	
Alondra-SS		Spatula	
Alondra-SD		Curved Dissector	
Alondra-AL	Trumpet 8° Angled	L Hook	
Alondra-AJ		J Hook	
Alondra-AS		Spatula	
Alondra-AD		Curved Dissector	



Company Information

Meta Biomed Co.,Ltd.

270 Osongsaengmyeong1-ro,
Osong-eup, Heungdeok-gu,
Cheongju-si, Chungbuk, Korea

- Phone: +82 43 230 8846
- Fax: + 82 43 230 8969
- boram@meta-biomed.com

www.meta-biomed.com

Meta Biomed Inc.

121 Independence Lane
Chalfont, PA 18914, USA

- Phone: +1 267 282 5893
- america@metabiomed-inc.com

www.buymetabiomedonline.com

Meta Biomed Europe

Haus der Wirtschaft, Raum 3.11
Wiesenstraße 35
45473 Mülheim an der Ruhr, Germany

- Phone: +49 208 309 9190
- europe@metabiomed-inc.com

www.buymetabiomedonline.com

H0-22-0